

Elements of a Preventative Maintenance Servicing

Meets or exceeds all College of American Pathology (CAP)TM preventative maintenance service requirements. (“WORDTM” version downloadable at MicroscopyUSA.com.)

- Evaluation**
Complete a performance check of all optical and mechanical functions; notify the laboratory staff of major concerns prior to servicing.
- Tear Down and Cleaning**
Carefully remove the oculars, head assembly, objectives (maintain sequence), slide holder, mechanical stage (if practical), condenser, blue / interference / neutral density or polarizing filters, field opening bracket with iris (if removable), and other accessories. Clean all optical surfaces using proper technique. Clean all non-optical solid surfaces including the body of the microscope using proper technique.
- Lubrication**
If needed, apply appropriate lubrication to all moving parts ... mechanical stage bearing track and rack and pinion gears, slide holder, focus mechanism, rotating nosepiece, condenser rack and pinion gears.
- Re-assemble**
- Kohler Alignment**
Complete Kohler Alignment using settings for the condenser alignment screws and the condenser and field iris diaphragms.
- Centration Error Check**
Check and record the specific values for Centration Error ... both horizontal and vertical.
- Calibrate Ocular Micrometer – only if requested by the clinical laboratory** – Use NIST* Stage Micrometer (*National Institute of Standardization of Technology)
- Full Function Check**
 - Parfocality, Parcentration and Resolution of all Objectives
 - Confirm that all ocular surfaces are clean by rotating oculars under high magnification.
 - Check the function of the Inter-Pupillary Distance (IPD) mechanism.
 - Check the function of focusing oculars and / or focusing eye-tubes.
 - Confirm stage collimation by focus stability during slide movement.
 - Check mechanical stage motion and control knob tension adjustments.
 - Check bulb socket alignment and condition.
 - Confirm proper focusing.
 - Check all accessories for proper function and alignment: phase contrast alignment, polarization and red compensator function, UV bulb alignment, micrometer calibration and documentation, dual-view system alignments, pointer bulbs, camera / CCD function.
 - Electrical Safety Check
- Calibrate Ocular Micrometer if Requested by Laboratory Staff**
- Notification of Concerns to Staff / Written Documentation to Lab Manager ...** conforming to CAP inspection requirements
- End-of-Day Verification** – a walk around of all lab departments to confirm that all microscopes are working well and that any adjustments are acceptable. **Total Time: ~ 35 minutes with experience.**